

Platelet Incubator (Room) Temperature Log



Month & Year: _____

If the temperature is not between +20 to +24°C please notify _____ immediately

If a unit is within 2 days of expiry, please notify _____ stock control

Date	Time AM	Time PM	Min T° +20°C	Max T° +24°C	Units agitating?	Thermometer reset?	Note (discolouration or visible contamination)
1 st							
2 nd							
3 rd							
4 th							
5 th							
6 th							
7 th							
8 th							
9 th							
10 th							
11 th							
12 th							
13 th							
14 th							
15 th							
16 th							
17 th							
18 th							
19 th							
20 th							
21 st							
22 nd							
23 rd							
24 th							
25 th							
26 th							
27 th							
28 th							
29 th							
30 th							
31 st							

FRM/SIS/67/01 – Platelet units should be stored lying on a platelet agitator in a temperature-controlled room or platelet incubator with a minimum/maximum thermometer placed close to the units. The temperature should remain between +20 and +24°C. The units should be continuously agitated and checked for discolouration or the presence of particles that may render them unusable.

Platelet Incubator (Room) Temperature Log



Month & Year: _____

If the temperature is not between +20 to +24°C please notify _____ immediately

If a unit is within 2 days of expiry, please notify _____ stock control

Date	Time AM	Time PM	Min T° +20°C	Max T° +24°C	Units agitating?	Thermometer reset?	Note (discolouration or visible contamination)
1 st							
2 nd							
3 rd							
4 th							
5 th							
6 th							
7 th							
8 th							
9 th							
10 th							
11 th							
12 th							
13 th							
14 th							
15 th							
16 th							
17 th							
18 th							
19 th							
20 th							
21 st							
22 nd							
23 rd							
24 th							
25 th							
26 th							
27 th							
28 th							
29 th							
30 th							
31 st							

FRM/SIS/67/01 – Platelet units should be stored lying on a platelet agitator in a temperature-controlled room or platelet incubator with a minimum/maximum thermometer placed close to the units. The temperature should remain between +20 and +24°C. The units should be continuously agitated and checked for discolouration or the presence of particles that may render them unusable.